



# BARON<sup>®</sup>

Platelet Rich Plasma (PRP),  
Naturally yours



Autologous Platelet Rich Plasma System

# BARON<sup>®</sup>



Produced for :  
AmericanCreations  
401 52nd St.  
West Des Moines , IA , USA  
Info@american-creations.com

☎ 25713133  
✉ info@pharmacliqueintl.com

THE VITALITY COMES FROM YOUR OWN CELLS

# BARON-Kit® - Product Highlights

## Simplicity:

Only a few steps are required to prepare 4 to 5 ml of Platelet-Rich Plasma (PRP).

## Low Blood Volume Required:

BARON-Kit® is designed to prepare Platelet Rich Plasma from only 8 ml of blood.

## Optimum Platelet Recovery:

BARON-Kit® you recover around 84% of the platelets from the whole blood. Platelets are concentrated in Autologous Plasma.

## Consistent Isolation of PRP:

BARON-Kit® uses a Non sticky physical barrier (thixotropic gel) to separate the red blood cells from the plasma and the cellular sediment. The preparation is safe and consistent With perfect controlled gel density.

## Red Blood Cell Depletion:

The level of Red Blood Cells is drastically reduced, with more than 99,7% of complete Red Blood Cell removal.



## Application:

Skin treatment, beauty industry, hair loss, osteoarthritis.

## Additive:

ACD/Sodium citrate+speration gel

## Advantages:

Crystal Tubes  
Sterile and pyrogen free  
Maximum active playlets with original morphology

## PRP volume per tube:

4ml-5ml

## Centrifusion process:

(3200) RPM for (6) minutes

BARON-Kit® contains (4) separate sterilized tubes

## Shelf life: 2 years

Blood Vol Per tube	PRP Vol Per tube	Platelet Recovery	Red blood cell Remanant	Platelet concentration factor
8ml	4 to 5ml	84%	<0.3%	10 x

